



Betamix AdhesionBarrierGel

Hyaluronic Acid
Prevents Adhesion
Fully Absorbable



Betamix AdhesionBarrierGel

Betamix is a clear, single use flowable gel that's used during different surgeries to help prevent adhesions, or internal scar tissue. Betamix acts as a temporary, protective barrier to separate tissues and reduce fibrosis and the formation of post-surgical adhesions.

Betamix is made up of Hyaluronic Acid which occurs naturally in the body and is used in pharmaceuticals, foods, and cosmetics. Once the gel is placed inside the body, it protects against adhesions for two weeks.

Betamix is non-pyrogenic and contains no animal components. No color additives are used in gel.



Advantages

- Ready to use with precision applicator
- Thorough coverage
- Colorless
- Prevents or reduces adhesion formation
- Safe, natural and absorbable
- Quick and easy application
- Cost effective

Ordering Information

Code	Volume
ABG01	1 ml
ABG105	1,5 ml
ABG02	2 ml
ABG03	3 ml
ABG04	4 ml
ABG05	5 ml
ABG705	7,5 ml
ABG10	10 ml
ABG15	15 ml
ABG20	20 ml

Betamix Spine
Betamix Tendon&Nerve
Betamix AbdominoPelvic

Betatech Medikal Cihazlar
San. Müm. İç ve Dış
Tic. Ltd.Şti.

İkitelli Org. San. Bölgesi
Atatürk Oto San. Sitesi
22. Sk. Ünal İş Merkezi No: 9,
Başakşehir-İstanbul, TÜRKİYE

T. +90 212 671 0683
F. +90 212 671 0684
E. info@betatechmedical.com

BETATECH
MEDICAL.COM

