





BETAMIX VAGINAL SLING

length: 450 mm (+/- 5 mm), width: 10 mm (+/- 0,5 mm), thickness: 0,34 mm (+/- 0,05 mm). Knit mesh size: 1 x 1,25 mm. Thread stitch length: 145 mm (+/- 5 mm). Concentration of knit holes: 55% (+/- 5 %). Density: 45 g/m² (+/- 5 %).



BETAMIX POP4

(four-arm) dimensions (trapezoid): height: 80 mm (+/- 1 mm), bases 60 and 90 mm (+/- 1 mm). Arms length: 180 and 195 mm (+/- 5 mm), width: 10 mm (+/- 0,5 mm). Thickness: 0,34 mm (+/- 0,05 mm). Knit mesh size: 1 x 1,25 mm. Thread stitch length: 145 mm (+/- 5 mm). Concentration of knit holes: 55% (+/- 5 %). Density: 45 g/m² (+/- 5 %).



BETAMIX POP2

(two-arm) dimensions: height: 103 mm (+/- 1 mm), width: 63 mm (+/- 1 mm). Arms: length: 190 mm (+/- 1 mm), width: 10 mm (+/- 0,5 mm). Thickness: 0,34 mm (+/- 0,05 mm). Knit mesh size: 1 x 1,25 mm. Thread stitch length: 145 mm (+/- 5 mm). Concentration of knit holes: 55% (+/- 5 %). Density: 45 g/m² (+/- 5 %).